

# ONYXA

---

MEDICAL

*IMPORTANT NOTICE: This document represents a suggested supply framework only. Final selection, preparation, and use of all medications and supplies are determined solely by the treating physician based on clinical judgment, patient needs, facility policies, and applicable local, state, and federal regulations.*

## TACTiQ Procedure Support Kit

Suggested Supplies for Tumescant Anesthesia, Aseptic Preparation, Procedure, and Post Care

### I. TACTiQ-T | Tumescant Anesthesia Components

- Xylocaine 2% with 1:100,000 epinephrine
- Lidocaine 1%
- Epinephrine (Adrenaline) 1 mg/mL
- Tranexamic Acid (TXA) 100 mg/mL
- Sodium Bicarbonate 8.4%
- Bacteriostatic 0.9% Sodium Chloride
- 18G needles (any length)
- 21G needles (any length)
- 30G needle
- 22G x 2 3/4 inch blunt tip cannula
- 22G x 4 inch blunt tip cannula (body applications)
- 20 cc syringe
- 10 cc syringe
- 1 cc syringe

### II. TACTiQ-P | Patient Prep and Aseptic Field Setup

- Hibiclens
- Puracyn Hypochlorous Acid
- Sterile water
- Sterile bowl
- Sterile gloves
- Sterile OR towels
- Tray cover
- Patient drape
- Purple surgical marker
- Surgical marker remover

### III. TACTiQ-PR | Procedure Day Clinical Supplies

- Sterile 2x2 or 4x4 gauze
- Additional sterile gauze as needed per case

### IV. TACTiQ-R | Post Procedure Care and Support

- Toppe foam
- Compression wrap or elastic net dressing

Suggested Supply Framework | For In-Service Preparation Use Only